

Focus™, FreeStyle™, & FreeStyle™ 5

Portable Oxygen Concentrators with UltraSense™

ALARM SHEET

Focus		
Alarm Type	Alarm Signal(s)	Action
Low Battery Alarm	Single ½ second blinking amber light; single ½ second audible alarm, followed by a 5-second pause.	Focus batteries should be plugged in an AC/DC power outlet as the battery needs to be charged.
Battery Shut Down	Two ½ second blinking amber lights; two ½ second audible alarms, followed by a 5-second pause.	Focus batteries should be plugged in an AC/DC power supply as the battery needs to be charged.
Overdraw Alarm	Three ½ second blinking amber lights; three ½ second audible alarms, followed by a 5-second pause.	Patients respiratory rate exceeds capability of Focus. Very uncommon alarm.
No Inspiration/Cannula Disconnect (15 minute delay)	Continuous amber light, constant audible alarm.	Patient forgot to turn Focus off, cannula is disconnected, patient could be a mouth breather.
Malfunction (Low O ₂ , High/Low Pressure)	Three ½ second blinking red lights; three ½-second audible alarms, followed by a 5-second pause.	Call for service.
FreeStyle		
Solid Amber Light	Solid amber light.	Service Required (5,000 Hours): arrange for service in the near future – no need for emergency pick-up
Low Battery Warning	Blinking 25% amber battery light (light illuminates on for ½ second every 5 seconds).	FreeStyle should be plugged in an AC/DC power outlet as the Battery needs to be charged.
Low Battery Alarm	Constant amber 25% (light illuminates on for ½ second then off for ½ second), rapid beeping audible alarm.	FreeStyle should be plugged in an AC/DC power outlet as the Battery needs to be charged.
Cannula Disconnect Alarm (15 minute delay)	Continuous red light, constant audible alarm.	Patient forgot to turn FreeStyle off, cannula is disconnected, patient could be a mouth breather.
Overdraw Alarm	Intermittent red light, intermittent audible alarm ½ second on, ½ second off.	Patients respiratory rate exceeds capability of Freestyle. Very uncommon alarm.
Malfunction (Low O ₂ , High/Low Pressure)	Constant red light, intermittent audible alarm, ½ second on, ½ second off.	POC maintenance required; equipment pick-up is necessary.
FreeStyle 5		
Low Battery Warning	Blinking 25% amber battery light, single intermittent (single beep) audible alarm.	FreeStyle 5 should be plugged into an AC/DC power outlet as the Battery needs to be charged.
Low Battery Alarm (Shut Down)	Blinking amber 25% light, double intermittent (beep, beep).	Battery Shut-down FreeStyle 5 should be plugged into an AC/DC power outlet as the Battery needs to be charged.
Cannula Disconnect Alarm (15 minute delay)	Continuous yellow light, constant audible alarm.	Patient forgot to turn FreeStyle off, cannula is disconnected, patient could be a mouth breather.
Overdraw Alarm	Yellow visual intermittent (blinking) light, rapid intermittent audible alarm (beep, beep, beep)	Patients respiratory rate exceeds capability of Freestyle. Very uncommon alarm. Patient should reduce current activity.
Malfunction (Low O ₂ , High/Low Pressure)	Constant red visual light, rapid intermittent audible alarm (beep, beep, beep).	POC maintenance required; equipment pick-up is necessary.

Contact Customer Support 1-800-482-2473 or Technical Service 770-721-7759